

Product Datasheet

Human SARS-CoV-2 (Covid-19) Qualitative ELISA Kit [Temporarily Out of Stock] orb638901

Description	Enzyme Immunoassay for the Qualitative Determination of SARS-CoV- 2 (Covid-19) Antigen in respiratory specimens and human serum.
Reactivity	Human
Sensitivity	refer IFU/Kit Insert
Tested Applications	ELISA
Target	SARS-CoV-2 (Covid-19) Qualitative
Storage	Store at 2-8°C
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	Precoated plates using human SARS-CoV-2 protein for capture/detection; type - sandwich ELISA