

Product Datasheet

RNA polII antibody orb640577

Description	Mouse monoclonal antibody to RNA polII
Species/Host	Mouse
Reactivity	Fungi, Human
Conjugation	Unconjugated
Tested Applications	ChIP, ELISA, FACS, IF, IHC-P, IP, WB
Immunogen	Peptide: 10 repeats of synthetic peptide YSPTSPS using chemically synthesized phospho-Ser5
Target	RNA polII
Note	For research use only.
Isotype	IgG1
Clonality	Monoclonal
Clone ID	CTD 4H8
Expiration Date	12 months from date of receipt.