

## Product Datasheet

### RSK1 p90 (Phospho-S380) antibody (orb650938)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | Rabbit polyclonal antibody to RSK1 p90 (Phospho-S380)...   |
| <b>Species/Host</b>        | Rabbit   |
| <b>Reactivity</b>          | Human, Mouse, Rat  |
| <b>Conjugation</b>         | Unconjugated   |
| <b>Tested Applications</b> | ICC, IF, IP, WB  |
| <b>Immunogen</b>           | Synthetic phospho-peptide corresponding to residues surrounding Ser380 of human RSK1.              |
| <b>Storage</b>             | Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. |
| <b>Note</b>                | For research use only.   |
| <b>MW</b>                  | 90 kDa   |
| <b>Uniprot ID</b>          | <a href="#">Q15418</a> , <a href="#">Q63531</a> , <a href="#">P18653</a>                           |
| <b>Dilution Range</b>      | WB:1:1, 000-1:5, 000ICC:1:50-1:200   |
| <b>Expiration Date</b>     | 12 months from date of receipt.  |