

Product Datasheet

Human proBNP (77-108) peptide orb81961

Description	pro-Brain Natriuretic Peptide, C-terminal (proBNP) (a.a. 77-108) is a purified cardiac marker.
Species/Host	Synthetic
Reactivity	Human
Conjugation	Unconjugated
Target	pro-Brain Natriuretic Peptide, C-terminal (proBNP) (a.a. 77-108). CXCR4, N terminal
Form/Appearance	Lyophilised.
Storage	Short to term store at -2 to -8°C. Long term store at to-20°C.Avoid multiple freeze/thaw cycles.
Note	For research use only.
Purity	Greater than 95% pure HPLC
Source	Synthetic
Biological Origin	Synthetic
Expiration Date	12 months from date of receipt.
Application Notes	Specific methodologies have not been tested using this product. Centrifuge before opening to ensure complete recovery of vial contents.