

Product Datasheet

CSFV E2 protein antibody {CF594} orb182960-CF594

Description	Rabbit polyclonal antibody to CSFV E2 protein
Species/Host	Rabbit
Reactivity	Virus
Conjugation	CF594
Tested Applications	ICC, IF
Immunogen	Full length CSFV E2 protein Recombinded
Target	CSFV E2 protein
Form/Appearance	Liquid: PBS, 2 mM Sodium Azide
Concentration	100 µg (in 411 µl): 0.24 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purified by Protein A
Dilution Range	To be determined by the end user
Expiration Date	12 months from date of receipt.