

Product Datasheet

E.coli MRJP1 protein orb244784

Description	Recombinant Apis mellifera Major royal jelly protein 1
Conjugation	Unconjugated
Target	MRJP1
Form/Appearance	Tris-based buffer, 50% glycerol
Storage	Store at -20°C. for extended storage, conserve at -20°C or -80°C.
Note	For research use only.
Protein Sequence	NILRGESLNKSLPILHEWKFFDYDFGSDERRQDAILSGEYDYKNNYPSDI DQWHDKIFVTMLRYNGVPSSLNVISKKVGDDGPLLQYPDWSFAKYDDCS GIVSASKLAIDKCDRLWVLDSGLVNNTQPMCSPKLLTFDLTTSQLLKQVE IPHDVAVNATTGKGRSSLAVQSLDCNTNSDTMVYIADEKGEGLIVYHNS DDSFHRLTSNTFDYDPKFTKMTIDGESYTAQDGISGMALSPMTNNLYYSP VASTS
Purity	>90% as determined by SDS-PAGE.
MW	74.2 kDa
Source	E. coli
Uniprot ID	O18330
Expiration Date	6 months from date of receipt.
Application Notes	This is GST-tag protein