

Product Datasheet

E.coli cry1Fb protein orb244790

Description	Recombinant Pesticidal crystal protein cry1Fb
Conjugation	Unconjugated
Target	cry1Fb
Form/Appearance	Tris-based buffer, 50% glycerol
Storage	Store at -20°C. for extended storage, conserve at -20°C or -80°C.
Note	For research use only.
Protein Sequence	VKGHVDVEEQNNHRSVLVPEWEAEVVSQEVVCPGRGYILRVTAYKEGYG EGCVTIHEVDNNTDELKFSNCEKEQVYPGNTVACNDYKNHGANACSSR NGGYDESYESNSSIPADYAPVYEEEEAYTDGQRGNPCEFNRGHTPLPAGYV TAELEYFPETDTVWVEIGETEGTFIVDSVELLLMEE
Purity	>90% as determined by SDS-PAGE.
MW	48.1 kDa
Source	E. coli
Uniprot ID	O66377
Expiration Date	6 months from date of receipt.
Application Notes	This is GST-tag protein