

Product Datasheet

Human ICT1 protein orb244801

Description	Recombinant human Peptidyl-tRNA hydrolase ICT1, mitochondrial
Species/Host	Human
Conjugation	Unconjugated
Target	ICT1
Form/Appearance	Tris-based buffer, 50% glycerol
Storage	Store at -20°C. for extended storage, conserve at -20°C or -80°C.
Note	For research use only.
Protein Sequence	LHKQKDGTEFKSIYSLDKLYPESQGSDTAWRVPNGAKQADSDIPLDRLTI SYCRSSGPGGQNVNKNVNSKAEVRFHLATAEWIAEPVRQKIAITHKNKINR LGELILTSESSRYQFRNLADCLQKIRDMITEASQTPKEPTKEDVKLHRIR IENMNRERLRQKRIHSAVKTSRRVDMD
Purity	>90% as determined by SDS-PAGE.
MW	36.4 kDa
Source	E. coli
Uniprot ID	Q14197
Expiration Date	6 months from date of receipt.
Application Notes	This is His-SUMO-tag protein