

Product Datasheet

Human JOSD1 protein orb244830

Description	Recombinant human Josephin-1
Species/Host	Human
Conjugation	Unconjugated
Target	JOSD1
Form/Appearance	Tris-based buffer, 50% glycerol
Storage	Store at -20°C. for extended storage, conserve at -20°C or -80°C.
Note	For research use only.
Protein Sequence	MSCVPWKGDKAKSESLELPQAAPPQIYHEKQRRELCALHALNNVFQDSNA FTRDTLQEIFQRLSPNTMVTPhkksMLGNGNYDVNVIMAALQTKGYEAVW WDKRRDVGVIAltnvmGfImnlPSSLCWGpLkLpLkRQHWICVREVGgAY YNLDSKlKMPeWIGGESElRkFlKHHLRGKNCeLLLVPEEVEAHQSWRT DV
Purity	>90% as determined by SDS-PAGE.
MW	39.2 kDa
Source	E. coli
Uniprot ID	Q15040
Expiration Date	6 months from date of receipt.
Application Notes	This is His-SUMO-tag protein