

Product Datasheet

Rat Hpx protein orb244277

Description	Recombinant rat Hemopexin
Species/Host	Rat
Conjugation	Unconjugated
Target	Hpx
Form/Appearance	Tris-based buffer, 50% glycerol
Storage	Store at -20°C. for extended storage, conserve at -20°C or -80°C.
Note	For research use only.
Protein Sequence	NPLPAAHETVAKGENGTKPDSVDVIEHCSDAWSFDATTMDHNGTMLFFKGE FVWRGHSGIRELISERWKNPVTSDAAFRGPDSVFLIKEDKVVVYPPEKK ENGYPKLFQEESPGIPYPPDAAVECHRGEQCSEGVLFFQGNRKWFWD FATRTQKERSWPAVGNCTAALRWLERYCYFQGNKFLRFNPVTGEVPPRYPLDA RDYFISCPGRGHGKLRNGTAHGNSTHPMHSRCNADPGLSALLSDHRGATY AFSGS
Purity	>90% as determined by SDS-PAGE.
MW	53 kDa
Source	E. coli
Uniprot ID	P20059
Expiration Date	6 months from date of receipt.
Application Notes	Full length of His-tag and expression region is 24-460aa