

Product Datasheet

Human ID2 protein orb244284

Description	Recombinant human DNA-binding protein inhibitor ID-2
Species/Host	Human
Conjugation	Unconjugated
Target	ID2
Form/Appearance	Tris-based buffer, 50% glycerol
Storage	Store at -20°C. for extended storage, conserve at -20°C or -80°C.
Note	For research use only.
Protein Sequence	MKAFSPVRSVRKNSLSDHSLGISRSKTPVDDPMSLLYNMND CYSK LKELV PSIPQNKVKMEILQHVIDYILDQLALDSHPTIVSLHHQRPGQNQASR TPLTTLNTDISILSLQASEFPSELMSNDSKALCG
Purity	>90% as determined by SDS-PAGE.
MW	30.9 kDa
Source	E. coli
Uniprot ID	Q02363
Expiration Date	6 months from date of receipt.
Application Notes	Full length of His-SUMO-tag and expression region is 1-134aa