

Product Datasheet

Mouse Hadh protein orb244848

Description	Recombinant mouse Hydroxyacyl-coenzyme A dehydrogenase, mitochondrial
Species/Host	Mouse
Conjugation	Unconjugated
Target	Hadh
Form/Appearance	Tris-based buffer, 50% glycerol
Storage	Store at -20°C. for extended storage, conserve at -20°C or -80°C.
Note	For research use only.
Protein Sequence	SSSSASAAAKKILIKHVTVIGGGLMGAGIAQVAAATGHTVVLDQTEDI LAKSKKGIIEESLKRMAKKKFTENPKAGDEFVEKTLSCSLSTSTDAASVVHS TDLVVEAIVENLKLKNELFQRDKFAAEHTIFASNTSSLQITNIANATTR QDRFAGLHFFNPVPMKLVKIKTPMITSQKTFESLVDFCKTLGKHPVSKK DTPGFIVNRLLPYLIEAVRLHERGDASKEDIDTAMKLGAGYPMGPFELL DYVGL
Purity	>90% as determined by SDS-PAGE.
MW	49 kDa
Source	E. coli
Uniprot ID	Q61425
Expiration Date	6 months from date of receipt.
Application Notes	This is His-SUMO-tag protein