

Product Datasheet

Rat Lgals3bp protein orb244340

Description	Recombinant rat Galectin-3-binding protein
Species/Host	Rat
Conjugation	Unconjugated
Target	Lgals3bp
Form/Appearance	Tris-based buffer, 50% glycerol
Storage	Store at -20°C. for extended storage, conserve at -20°C or -80°C.
Note	For research use only.
Protein Sequence	AKDGDMLRVNGASASEGRVEIFYRGRWGTVCDNLWNLDAHVVCRALGYE NATQALSRAAFGPGKGPIMLDEVECTGNESSLANCSSLGWMVSHCGHEKD AGVVCSNDSRGIHILDLSGELPDALGQIFDSQQDCDLFIQVTGQGHGDLS LCAHTLILRTNPEAQALWQVVGSSVIMRVDAECMPVVRDFLRYFYSSRIE VSMSSVKCLHKLASAYGATELQGYCGRLFVTLQPDPFHTPLELYEYQAQ ATGDS
Purity	>90% as determined by SDS-PAGE.
MW	65.9 kDa
Source	E. coli
Uniprot ID	070513
Expiration Date	6 months from date of receipt.
Application Notes	Full length of His-tag and expression region is 19-574aa