

Product Datasheet

Human TMX3 protein orb244914

Description	Recombinant human Protein disulfide-isomerase TMX3
Species/Host	Human
Conjugation	Unconjugated
Target	TMX3
Form/Appearance	Tris-based buffer, 50% glycerol
Storage	Store at -20°C. for extended storage, conserve at -20°C or -80°C.
Note	For research use only.
Protein Sequence	KGFVEDLDESEFKENRNDDIWLVDYAPWCGHCKKLESIWNEVGLEMKSIG SPVKVKGMDATSYSSIASEFGVRYPTIKLLKGDLAYNYRGPRTKDDIIE FAHRVSGALIRPLPSQQMFEHMQKRHRVFFVYVGGESPLKEYIDAASEL IVYTYFFSASEEVVPEVIFKI
Purity	>90% as determined by SDS-PAGE.
MW	35.6 kDa
Source	E. coli
Uniprot ID	Q96JJ7
Expiration Date	6 months from date of receipt.
Application Notes	This is His-SUMO-tag protein