

Product Datasheet

Human MYLIP protein orb244926

Description	Recombinant human E3 ubiquitin-protein ligase MYLIP
Species/Host	Human
Conjugation	Unconjugated
Target	MYLIP
Form/Appearance	Tris-based buffer, 50% glycerol
Storage	Store at -20°C. for extended storage, conserve at -20°C or -80°C.
Note	For research use only.
Protein Sequence	MLCYVTRPDAVLMEVEVEAKANGEDCLNQVCRRLGIIIVDYFGLQFTGSK GESLWLNLRNRISQQMDGLAPYRLKLRVKFFVPHLILQEQRHIFFLHI KEALLAGHLLCSPEQAVELSALLAQTKFGDYNQNTAKYNYEELCAKELSS ATLNSIVAKHKELEGTSQASAEYQVLQIVSAMENYGIWHSVRDSEGQKL LIGVGPEGISICKDDFSPINRIAYPVVQMATQSGKNVYLTVTKESGNSIV LLFKM
Purity	>90% as determined by SDS-PAGE.
MW	77.3 kDa
Source	E. coli
Uniprot ID	Q8WY64
Expiration Date	6 months from date of receipt.
Application Notes	This is GST-tag protein