

Product Datasheet

Rabbit Mfi2 protein orb244374

Description	Recombinant Rabbit Melanotransferrin
Species/Host	Rabbit
Conjugation	Unconjugated
Target	Mfi2
Form/Appearance	Tris-based buffer, 50% glycerol
Storage	Store at -20°C. for extended storage, conserve at -20°C or -80°C.
Note	For research use only.
Protein Sequence	SVEVRWCTASEPEQQKCEDMSQAFREAGLQPALLCVQGTSAADHCVQLIAA HEADAITLDGGAIYEAGKEHGLKPVVGEVYDQEVGTSYYAVAVVKKRSSNV TINTLRGVKSCHTGINRTVGNVNPVGYLVDSGRLSVMGCDVLKAVSEYFG GSCVPGAGETRYSESLCRLCRGDTSGEGVCDKSPLERYDYSGAFRCLAE GAGDVAFVKHSTVLENTDGRITLPSWGHMLMSRDFELLCRDGSRASVTEWQ HCHLA
Purity	>90% as determined by SDS-PAGE.
MW	91.6 kDa
Source	E. coli
Uniprot ID	O97490
Expiration Date	6 months from date of receipt.
Application Notes	Partial of Melanotransferrin of His-SUMO-tag and expression region is 20-710aa