

Product Datasheet

Rat Mmp7 protein orb244379

| | |
|--------------------------|--|
| Description | Recombinant rat Matrilysin |
| Species/Host | Rat |
| Conjugation | Unconjugated |
| Target | Mmp7 |
| Form/Appearance | Tris-based buffer, 50% glycerol |
| Storage | Store at -20°C. for extended storage, conserve at -20°C or -80°C. |
| Note | For research use only. |
| Protein Sequence | FSLMPNSPKWHSRTVTYRIVSYTTDLPRFLVDQIVKRALRMWSMQIPLNF KRVSWGTAIIIIGFARGDHGDNFPFDGPGNTLGHAFAPGPGGLGGDAHFDK DEYWTDGEDSGVNFVATHELGHSLGLGHSSVPSSVMYPTYQGDHSEDF SLTKDDIAGIQKLYGKRNL |
| Purity | >90% as determined by SDS-PAGE. |
| MW | 46.3 kDa |
| Source | E. coli |
| Uniprot ID | P50280 |
| Expiration Date | 6 months from date of receipt. |
| Application Notes | Full length of GST-tag and expression region is 98-267aa |