

Product Datasheet

Human DBNDD2 protein orb244945

Description	Recombinant human Dysbindin domain-containing protein 2
Species/Host	Human
Conjugation	Unconjugated
Target	DBNDD2
Form/Appearance	Tris-based buffer, 50% glycerol
Storage	Store at -20°C. for extended storage, conserve at -20°C or -80°C.
Note	For research use only.
Protein Sequence	MDPNPRAALERQQLRLRERQKFFEDILQPETEFVFPPLSHLHLESQRPPIG SISSMEVNVDTLEQVELIDLGDPAADVFLPCEDPPPTPQSSGVDNHLEE LSLPVPTSDRTTSRTSSSSSSDSSTNLHSPNPSDDGADTPLAQSDEEEER GDGGAEPGACS
Purity	>90% as determined by SDS-PAGE.
MW	33.48 kDa
Source	E. coli
Uniprot ID	Q9BQY9
Expiration Date	6 months from date of receipt.
Application Notes	This is His-SUMO-tag protein