

Product Datasheet

Human DKK2 protein orb244953

Description	Recombinant human Dickkopf-related protein 2
Species/Host	Human
Conjugation	Unconjugated
Target	DKK2
Form/Appearance	Tris-based buffer, 50% glycerol
Storage	Store at -20°C. for extended storage, conserve at -20°C or -80°C.
Note	For research use only.
Protein Sequence	KLNSIKSSLGGETPGQAANRSAGMYQGLAFGGSKKGKNLGQAYPCSSDKE CEVGRYCHSPHQSSACMVCRRKKKRCHRDGMCCPSTRCNGICIPVTES ILTPHIPALDGTQHRDRNHGHYSNHDLGWQNLGRPHTKMSHIKGGHEGDPC LRSSDCIEGFCCARHFWTKICKPVLHQGEVCTKQRKKGSHGLEIFQRDCD AKGLSCKVWKDATYSSKARLHVCQKI
Purity	>90% as determined by SDS-PAGE.
MW	41 kDa
Source	E. coli
Uniprot ID	Q9UBU2
Expiration Date	6 months from date of receipt.
Application Notes	This is His-SUMO-tag protein