

Product Datasheet

Human MYEF2 protein orb244394

Description	Recombinant human Myelin expression factor 2
Species/Host	Human
Conjugation	Unconjugated
Target	MYEF2
Form/Appearance	Tris-based buffer, 50% glycerol
Storage	Store at -20°C. for extended storage, conserve at -20°C or -80°C.
Note	For research use only.
Protein Sequence	GGMNRIGGGIGFGGLEAMNSMGGFGGVGRMGELYRGAMTSSMERDFGRGD IGINRFGGDSFGRLLGGGMGSMNSVTGGMGMGLDRMSSSFDRMGPGIGAIL ERSIDMDRGFLSGPMGSGMRERIGSKGNQIFVRNLPFDLTWQKLKEKFSQ CGHVMFAEIKMENGKSKGCGTVRFDSPESA EKACRIMNGIKISGREIDVR LDRNA
Purity	>90% as determined by SDS-PAGE.
MW	37.8 kDa
Source	E. coli
Uniprot ID	Q9P2K5
Expiration Date	6 months from date of receipt.
Application Notes	Partial of His-SUMO-tag and expression region is 372-576aa