

Product Datasheet

E.coli At2g24580 protein [Out of Stock] orb244959

Description	Recombinant Arabidopsis thaliana Probable sarcosine oxidase
Conjugation	Unconjugated
Target	At2g24580
Form/Appearance	Tris-based buffer, 50% glycerol
Storage	Store at -20°C. for extended storage, conserve at -20°C or -80°C.
Note	For research use only.
Protein Sequence	MEYSDDGRFDVIVVGAGVMGSSAAYQLAKRGQKTLLEQFDLHHRGSSH GESRTIRATYPEDYYSMVSESTRLWAAAQSEIGYKVHFPTQQFDMGPAD QQSLLSVVATCQKHGLAHRVMDSHAVSEHFSGRISIPENWIGVSTELGGI IKPTKAVSMFQTLAIGHGAILRDNTKVANIKRDGESGEGVIVCTVKGDKF YGKKCIVTAGAWISKLVKTVAGIDFPVEPLETTVCYWRIKEGHEEKFTID GEFPT
Purity	>90% as determined by SDS-PAGE.
MW	49.8 kDa
Source	E. coli
Uniprot ID	Q9SJA7
Expiration Date	6 months from date of receipt.
Application Notes	This is His-tag protein