

Product Datasheet

Human MAF1 protein orb244962

| | |
|--------------------------|---|
| Description | Recombinant human Repressor of RNA polymerase III transcription MAF1 homolog |
| Species/Host | Human |
| Conjugation | Unconjugated |
| Target | MAF1 |
| Form/Appearance | Tris-based buffer, 50% glycerol |
| Storage | Store at -20°C. for extended storage, conserve at -20°C or -80°C. |
| Note | For research use only. |
| Protein Sequence | MKLENSFEAINSQLTVETGDAHIIGRIESYSCKMAGDDKHMFKQFCQE GQPHVLEALSPPQTSGLSPSRLSKSQGGEEGPLSDKSRKTLFYLIATL NESFRPDYDFSTARSHFEFSREPSLSWVVNAVNCSLFSAVREDFKDLKPQL WNAVDEEICLAECDIYSYNDLSDPFGEDGSLWSFNFFYNKRLKRIVF FSCRSISGSTYTPSEAGNELDMELGEEVEEESRSRSGAEETSTMEEDR VPVIC |
| Purity | >90% as determined by SDS-PAGE. |
| MW | 44.9 kDa |
| Source | E. coli |
| Uniprot ID | Q9H063 |
| Expiration Date | 6 months from date of receipt. |
| Application Notes | This is His-SUMO-tag protein |