

Product Datasheet

Human NOL3 protein orb244408

Description	Recombinant human Nucleolar protein 3
Species/Host	Human
Conjugation	Unconjugated
Target	NOL3
Form/Appearance	Tris-based buffer, 50% glycerol
Storage	Store at -20°C. for extended storage, conserve at -20°C or -80°C.
Note	For research use only.
Protein Sequence	MGNAQERPSETIDRERKRLVETLQADSGLLLDALLARGVLTGPEYEALDA LPDAERRVRRLLLLVQKGGEAACQELLRCAQRTAGAPDPAWDWQHVGPGY RDRSYDPPCPGHWTPPEAPGSGTTCPGLPRASDPDEAGGPEGSEAVQSGTP EEPEPELEAEASKEAEPEPEPEPELEPEAEAEPEPELEPEPDPEPEPDFE ERDESEDS
Purity	>90% as determined by SDS-PAGE.
MW	38.6 kDa
Source	E. coli
Uniprot ID	060936
Expiration Date	6 months from date of receipt.
Application Notes	Full length of Isoform 2 of His-SUMO-tag and expression region is 1-208aa