

Product Datasheet

Human OLA1 protein orb244421

Description	Recombinant human Obg-like ATPase 1
Species/Host	Human
Conjugation	Unconjugated
Target	OLA1
Form/Appearance	Tris-based buffer, 50% glycerol
Storage	Store at -20°C. for extended storage, conserve at -20°C or -80°C.
Note	For research use only.
Protein Sequence	MPPKKGDDGIKPPPIIGRFGTSLKIGIVGLPNVGKSTFFNVLNLSQASAE NFPFCTIDPNESRVPVDERFDLCQYHKPASKIPAFNLNVVDIAGLVKGA HNGQGLGNAFLSHISACDGIHLTRAFEDDDITHVEGSVDPIRDIEIHE ELQKDEEMIGPIIDKLEKVAVRGGDKKLPKPEYDIMCKVKSVIDQKKPV RFYHDWWDKEIEVLNKHFLTSKPMVYLVNLSEKDYIRKKNKWLIIKEW VDKYD
Purity	>90% as determined by SDS-PAGE.
MW	60.7 kDa
Source	E. coli
Uniprot ID	Q9NTK5
Expiration Date	6 months from date of receipt.
Application Notes	Full length of His-SUMO-tag and expression region is 1-396aa