

## Product Datasheet

### Mouse Pdia2 protein orb244433

<b>Description</b>	Recombinant mouse Protein disulfide-isomerase A2
<b>Species/Host</b>	Mouse
<b>Conjugation</b>	Unconjugated
<b>Target</b>	Pdia2
<b>Form/Appearance</b>	Tris-based buffer, 50% glycerol
<b>Storage</b>	Store at -20°C. for extended storage, conserve at -20°C or -80°C.
<b>Note</b>	For research use only.
<b>Protein Sequence</b>	QGEEPGGPSEVLPEEPTGEEVPKEDGILVLNHRTLSSLALQEHSALMVEFY APWCGHCKELAPEYSKAAALLAAESAVVTLAKVDGPAEPELTKEFEVVG PTLKFFQNGNRTNPEEYAGPKTAEGIAEWLRRRVGPSATHLEDEEGVQAL MAKWDMMVIGFFQDLQKDMATFLALAKDALDMTFGFTDQPQLFEKFGLT KDTVVLFFKKFDEGRADFPVDKETGLDLGDLRSRFLVIHSMHLVTEFNSQTS PKIFA
<b>Purity</b>	>90% as determined by SDS-PAGE.
<b>MW</b>	72.2 kDa
<b>Source</b>	E. coli
<b>Uniprot ID</b>	<b>D3Z6P0</b>
<b>Expiration Date</b>	6 months from date of receipt.
<b>Application Notes</b>	Full length of His-SUMO-tag and expression region is 21-527aa