

Product Datasheet

Human NAGK protein orb244981

Description	Recombinant human N-acetyl-D-glucosamine kinase
Species/Host	Human
Conjugation	Unconjugated
Target	NAGK
Form/Appearance	Tris-based buffer, 50% glycerol
Storage	Store at -20°C. for extended storage, conserve at -20°C or -80°C.
Note	For research use only.
Protein Sequence	AAIYGGVEGGGTRSEVLLVSEDGKILAEADGLSTNHWLIGTDKCOVERINE MVNRAKRKAGVDPLVPLRSLGLSLSGGDQEDAGRILIEELRDRFPYLSES YLITTTDAAGSIATATPDGGVVLISGTGSGNCRNLINPDGSESGCGGWHMMG DEGSAYWIAHQAVKIVFDSIDNLEAAPHDIGYVKQAMFHYFQVPDRLGIL THLYRDFDKCRFAGFCRKAIEGAQQGDPLSR YIFRKAGEMLGRHIVAVLP EIDPV
Purity	>90% as determined by SDS-PAGE.
MW	53.2 kDa
Source	E. coli
Uniprot ID	Q9UJ70
Expiration Date	6 months from date of receipt.
Application Notes	This is His-SUMO-tag protein