

Product Datasheet

Human CDC23 protein orb244991

Description	Recombinant human Cell division cycle protein 23 homolog
Species/Host	Human
Conjugation	Unconjugated
Target	CDC23
Form/Appearance	Tris-based buffer, 50% glycerol
Storage	Store at -20°C. for extended storage, conserve at -20°C or -80°C.
Note	For research use only.
Protein Sequence	LSINSDFSDLREIKKQLLLIAGLTRERGLLHSSKWSAELAFSLPALPLAE LQPPPIITEEDAQDMDAYTLAKAYFDVKEYDRAAHFLHGCNSKKAYFLYM YSRYLSGEEKKKDDETVDSLGPLEKGQVKNEALRELRVLSKKHQARELDG FGLYLYGVVLRKLDLVKEAIDVFVEATHVLPPLHWGAWLELCNLITDKEML KFLSLPDTWMKEFFLAHIYTELQLIEEALQKYQNLIDVGFSSSYIVSQI AVAYH
Purity	>90% as determined by SDS-PAGE.
MW	71.1 kDa
Source	E. coli
Uniprot ID	Q9UJX2
Expiration Date	6 months from date of receipt.
Application Notes	This is His-tag protein