

Product Datasheet

Human RELN protein orb244491

Description	Recombinant human Reelin
Species/Host	Human
Conjugation	Unconjugated
Target	RELN
Form/Appearance	Tris-based buffer, 50% glycerol
Storage	Store at -20°C. for extended storage, conserve at -20°C or -80°C.
Note	For research use only.
Protein Sequence	AAGYYPRFSPFFFLCTHHGELEGDGEQGEVLISLHIAGNPTYVPGQEYH VTISTSTFFDGLLVTGLYTSTSVQASQSIGGSSAFGFGIMSDHQFGNQFM CSVVASHVSHLPTTNLSFIWIAPPAGTGCVNFMATATHRGQVIFKDALAQ QLCEQGAPTDVTVHPLAEIHSDSILRDDFDSYHQLQLNPNIWVECNC ETGEQCGAIMHGNAVTFCEPYGPRELITT
Purity	>90% as determined by SDS-PAGE.
MW	29 kDa
Source	E. coli
Uniprot ID	P78509
Expiration Date	6 months from date of receipt.
Application Notes	Partial of the full length of 26-3460aa of His-tag and expression region is 26-254aa