

Product Datasheet

Robo antibody {CF647} orb100869-CF647

Description	Rabbit Polyclonal to Robo.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Predicted Reactivity	Rat
Conjugation	CF647
Tested Applications	ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from human Robo: 51-150/1651
Form/Appearance	Liquid
Concentration	100 µg in 500 µl (0.20 mg/ml)
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
MW	181kDa
Uniprot ID	Q9Y6N7, O89026
Dilution Range	To be determined by the end user
Expiration Date	12 months from date of receipt.