

Product Datasheet

Human SRSF1 protein orb244537

Description	Recombinant human Serine/arginine-rich splicing factor 1
Species/Host	Human
Conjugation	Unconjugated
Target	SRSF1
Form/Appearance	Tris-based buffer, 50% glycerol
Storage	Store at -20°C. for extended storage, conserve at -20°C or -80°C.
Note	For research use only.
Protein Sequence	SGGGVIRGPAGNNDCRIYVGNLPPDIRTKDIEDVFYKYGAIRDIDLKNRR GGPPFAFVEFEDPRDAEDAVYGRDGYDYDGYRLRVEFPRSGRGTGRGGGG GGGGGAPRGRYGPPSRSENRVVVSGLPPSGSWQDLKDHMREAGDVCYAD VYRDGTGVVEFVRKEDMTYAVRKLNTKFRSHEGETAYIRVKVDGPRSPS YGRSRSRSRSRSRSRSNSRSRSYSPRRSRGSPRYSRHSRHSRST
Purity	>90% as determined by SDS-PAGE.
MW	43.6 kDa
Source	E. coli
Uniprot ID	Q07955
Expiration Date	6 months from date of receipt.
Application Notes	Full length of His-SUMO-tag and expression region is 2-248aa