

Product Datasheet

Human DAXX protein orb245070

Description	Recombinant human Death domain-associated protein 6
Species/Host	Human
Conjugation	Unconjugated
Target	DAXX
Form/Appearance	Tris-based buffer, 50% glycerol
Storage	Store at -20°C. for extended storage, conserve at -20°C or -80°C.
Note	For research use only.
Protein Sequence	MATANSIIVLDDDDDEDEAAAQPGPSHLPNAASPGAEAPSSSEPHGARGS SSSGGKKCYKLENEKLFEEFLELCKMQTADHPEVVPFLYNRQQRHSLFL ASAEFCNILSRVLSRARSRAKLYVYINELCTVLKAHSAKKKLNLA PAAT TSNEPSGNNPPTHLSLDPTNAENTASQSPRTRGSRRQIRLEQLLALYVA EIRRLQEKELDLSELDDPDSAYLQEARLKRKLI
Purity	>90% as determined by SDS-PAGE.
MW	53.1 kDa
Source	E. coli
Uniprot ID	Q9UER7
Expiration Date	6 months from date of receipt.
Application Notes	This is GST-tag protein