

Product Datasheet

Human YWHAB protein orb245088

Description	Recombinant human 14-3-3 protein beta/alpha
Species/Host	Human
Conjugation	Unconjugated
Target	YWHAB
Form/Appearance	Tris-based buffer, 50% glycerol
Storage	Store at -20°C. for extended storage, conserve at -20°C or -80°C.
Note	For research use only.
Protein Sequence	MTMDKSELVQKAKLAEQAERYDDMAAAMKAVTEQGHLSNEERNLLSVAY KNVVGARRSSWRVISSIEQKTERNEKKQOMGKEYREKIEAELQDICNDVL ELLDKYLIPNATQPESKVFYLMKMGDYFRYLSEVASGDNKQTTVSNSQQA YQEAFEISKKEMQPTHPIRLGLALNFSVFYIEILNSPEKACSLAKTAFDE AIAELDTLN EESYKDSTLIMQLLRDNLTLWTSENQGDGEGEN
Purity	>90% as determined by SDS-PAGE.
MW	55.5 kDa
Source	E. coli
Uniprot ID	P31946
Expiration Date	6 months from date of receipt.
Application Notes	This is GST-tag protein