

Product Datasheet

Human SYT1 protein orb244571

Description	Recombinant human Synaptotagmin-1
Species/Host	Human
Conjugation	Unconjugated
Target	SYT1
Form/Appearance	Tris-based buffer, 50% glycerol
Storage	Store at -20°C. for extended storage, conserve at -20°C or -80°C.
Note	For research use only.
Protein Sequence	KNAINMKDVKDLGKTMKDQALKDDDAETGLTDGEEKEEPKEEEKLGKQLQY SLDYDFQNNQLLVGIIQAAELPALDMGGTSDPYVKVFLLPDKKKKFETKV HRKTLNPFVNEQFTFKVPYSELGGKTLVMAVYDFDRFSKHDIIGEFKVPM NTVDFGHVTEEWRDLQSAEKEEQEKLGDICFSLRYVPTAGKLTVVILEAK NLKKMDVGGGLSDPYVKIHLMQNGKRLKKKKTTIKNTLNPYYNESFSFEV PFEQI
Purity	>90% as determined by SDS-PAGE.
MW	40.4 kDa
Source	E. coli
Uniprot ID	P21579
Expiration Date	6 months from date of receipt.
Application Notes	Cytoplasmic domain of His-tag and expression region is 99-416aa