

Product Datasheet

Human C14orf3 protein orb245243

Description	Recombinant human Activator of 90 kDa heat shock protein ATPase homolog 1
Species/Host	Human
Conjugation	Unconjugated
Target	C14orf3
Form/Appearance	Tris-based buffer, 50% glycerol
Storage	Store at -20°C. for extended storage, conserve at -20°C or -80°C.
Note	For research use only.
Protein Sequence	MAKWGEGDPRWIVEERADATNVNNWHWTERDASNWSTDKCLKTLFLAVQVQ NEEGKCEVTEVSKLDGEASINNRKKGKLIFFYEWVSVKLNWTGTSKSGVQYK GHVEIPNLSDENSVDEVEISVSLAKDEPDTNLVALMKEEGVKLLREAMGI YISTLKTEFTQGMILPTMNGESVDPVGPALKTEERKAKPAPSKTQARPV GVKIPTCKITLKETFLTSPPELYRVFTTQELVQAFTHAPATLEADRGGKF HMVDG
Purity	>90% as determined by SDS-PAGE.
MW	65.7 kDa
Source	E. coli
Uniprot ID	O95433
Expiration Date	6 months from date of receipt.
Application Notes	This is GST-tag protein