

Product Datasheet

Human CPA3 protein orb245799

Description	Recombinant human Mast cell carboxypeptidase A
Species/Host	Human
Conjugation	Unconjugated
Target	CPA3
Form/Appearance	Tris-based buffer, 50% glycerol
Storage	Store at -20°C. for extended storage, conserve at -20°C or -80°C.
Note	For research use only.
Protein Sequence	IPGRHSYAKYNNWEKIVAWTEKMMDKYPEMVSRIKIGSTVEDNPLYVLKI GEKNERRKAIFTDCGIHAREWVSPAFCQWFVYQATKTYGRNKIMTKLLDR MNFYILPVFNVDGYIWSWTKNRMWRKNRSKNQNSKICIGTDLNRNFNASWN SIPNTNDPCADNYRGSAPSEKETKAVTNFIRSHLNEIKVYITFHSYSQM LLFPYGYTSKLPPNHEDLAKVAKIGTDVLS TRYETRYIYGPIESTIYPIS GSSLD
Purity	>90% as determined by SDS-PAGE.
MW	38 kDa
Source	Yeast
Uniprot ID	P15088
Expiration Date	6 months from date of receipt.
Application Notes	This is His-tag protein