

Product Datasheet

Human CSNK1E protein orb245802

Description	Recombinant human Casein kinase I isoform epsilon
Species/Host	Human
Conjugation	Unconjugated
Target	CSNK1E
Form/Appearance	Tris-based buffer, 50% glycerol
Storage	Store at -20°C. for extended storage, conserve at -20°C or -80°C.
Note	For research use only.
Protein Sequence	MELRVGNKYRLGRKIGSGSFGDIYLGANIASGEEVAIKLECVKTKHPQLH IESKFYKMMQGGVGIPSIKWCGAEGDYNVMVMELLGPSLEDLNFCSRKF SLKTVLLLADQMISRIEYIHSKNFIHRDVKPDNFLMGLGKKGNLVYIIDF GLAKKYRDARTHQHIPPYRENKNLTGTARYASINTHLGIEQSRDDLES LG YVLMYFNLGSLPWQGLKAATKRQKYERISEKKMSTPIEVLCGYPSEFST YLNFC
Purity	>90% as determined by SDS-PAGE.
MW	49.3 kDa
Source	Yeast
Uniprot ID	P49674
Expiration Date	6 months from date of receipt.
Application Notes	This is His-tag protein