

Product Datasheet

Human FCGRT protein [Out of Stock] orb245835

Description	Human FCGRT protein
Species/Host	Human
Conjugation	Unconjugated
Target	Human FCGRT
Form/Appearance	Tris-based buffer, 50% glycerol
Storage	Store at -20°C. for extended storage, conserve at -20°C or -80°C.
Note	For research use only.
Protein Sequence	AESHLSLLYHLTAVSSPAGTPAFWVSGWLGPPQYLSYNSLRGEAEP CGA WVWENQVSWYWEKETTDLRIKEKLFLEAFKALGGKGPYTLQGLL GCELGP DNTSVPTAKFALNGEEFMNFDLKQGTWGGDWPEALAI SORWQQDKAANK ELTFLFSCPHRLREHLERGRGNLEWKEPPSM RLKARPSSPGFSVLTCSA FSFYPPQLRFLRNLAAAGTGQGD FGPNSDGSFHASSSLTVKSGDEHHY CCIVQ
Purity	>90% as determined by SDS-PAGE.
MW	32.4 kDa
Source	Yeast
Uniprot ID	P55899
Expiration Date	6 months from date of receipt.
Application Notes	This is His-tag protein