

Product Datasheet

Mouse Gcgr protein orb245845

Description	Recombinant mouse Glucagon receptor
Species/Host	Mouse
Conjugation	Unconjugated
Target	Gcgr
Form/Appearance	Tris-based buffer, 50% glycerol
Storage	Store at -20°C. for extended storage, conserve at -20°C or -80°C.
Note	For research use only.
Protein Sequence	AQVMDFLFEKWKLYSDQCHHNLSELLPPPELVCNRTFDKYSCWPDTTPNT TANISCPWYLPWYHKVQHRLVFKRCGPDGQWVRGPRGQPWRNASQCQLDD EEIEVQKGVAKMYSSQQ
Purity	>90% as determined by SDS-PAGE.
MW	15.8 kDa
Source	Yeast
Uniprot ID	Q61606
Expiration Date	6 months from date of receipt.
Application Notes	This is His-tag protein