

Product Datasheet

Human HOXB4 protein orb245865

Description	Recombinant human Homeobox protein Hox-B4
Species/Host	Human
Conjugation	Unconjugated
Target	HOXB4
Form/Appearance	Tris-based buffer, 50% glycerol
Storage	Store at -20°C. for extended storage, conserve at -20°C or -80°C.
Note	For research use only.
Protein Sequence	MAMSSFLINSNYVDPKFPPEEYSQSDYLPDHPGYYAGGQRRESSFQP EAGFGRRRAACTVQRYAACRDPGPPPPPPPPPPPPPPGLSPRAPAPPPAG ALLPEPGQRCEAVSSPPPPCAQNPLHPSPSHSACEPVVYPWMRKVHV STVNPNYAGGEPKRSRTAYTRQQVLELEKEFHYNRYLTRRRRVEIAHALC LSERQIKIWFQNRMMKWKDKHKL PNTKIRSGGAAGSAGGPPGRPNGGPRA L
Purity	>90% as determined by SDS-PAGE.
MW	29.6 kDa
Source	Yeast
Uniprot ID	P17483
Expiration Date	6 months from date of receipt.
Application Notes	This is His-tag protein