

Product Datasheet

Human Thrombospondin-1 protein orb245326

Description	Recombinant human Thrombospondin-1 protein
Species/Host	Human
Conjugation	Unconjugated
Target	Thrombospondin-1
Form/Appearance	Tris-based buffer, 50% glycerol
Storage	Store at -20°C. for extended storage, conserve at -20°C or -80°C.
Note	For research use only.
Protein Sequence	NRIPESGGDNSVFDIFELTGAARKGSGRRLVKGPDPPSSPAFRIEDANLIP PVPDDKFQDLVDAVRAEKGFLLLASLRQMKKTRGTLALERKDHSGQVFS VVSNGKAGTLDLSLTVQKGQHVVSVEEALLATGQWKSITLFVQEDRAQLY IDCEKMENAELDVPIQSVFTRDLASIALRLIAKGGVNDNFQGVLQNVRFV FGTTPEDILRNKGCSSSTSVLLTLDNNVWNGSSPAIRTNYIGHKTKDLQA ICGIS
Purity	>90% as determined by SDS-PAGE.
MW	40 kDa
Source	E. coli
Uniprot ID	P07996
Expiration Date	6 months from date of receipt.
Application Notes	This is His-tag protein