

Product Datasheet

Human IDUA protein orb245872

Description	Recombinant human Alpha-L-iduronidase
Species/Host	Human
Conjugation	Unconjugated
Target	IDUA
Form/Appearance	Tris-based buffer, 50% glycerol
Storage	Store at -20°C. for extended storage, conserve at -20°C or -80°C.
Note	For research use only.
Protein Sequence	APHLVHVDAARALWPLRRFWRSTGFCPLPHSQADQYVLSWDQQLNLAYV GAVPHRGIKQVRTHWLELVTTTRGSTGRGLSYNFTHLDGYLDLLRENQLL PGFELMGSASGHFTDFEDKQQVFEWKDLVSSLARRYIGRYGLAHVSKWNF ETWNEPDHHDNDVSMQMGFLNYYDACSEGLRAASPALRLGGPGDSFHT PPRSPLSWGLLRHCHDGTNFFTGEAGVRLDYISLHRKRGARSSISILEQEK VVAQQ
Purity	>90% as determined by SDS-PAGE.
MW	71.9 kDa
Source	Yeast
Uniprot ID	P35475
Expiration Date	6 months from date of receipt.
Application Notes	This is His-tag protein