

Product Datasheet

H5N6 hemagglutinin antibody {CF647} orb221432-CF647

Description	Rabbit polyclonal antibody to HDGF
Species/Host	Rabbit
Conjugation	CF647
Tested Applications	ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from the region aa361-440 of human H5N6 Hemagglutinin
Target	H5N6 hemagglutinin
Form/Appearance	Liquid: PBS, 2 mM Sodium Azide
Concentration	100 µg (in 411 µl): 0.24 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purified by Protein A
MW	63 kDa
Dilution Range	To be determined by the end user
Expiration Date	12 months from date of receipt.