

## Product Datasheet

### SARS-CoV-2 Nucleocapsid protein orb669318

<b>Description</b>	Recombinant SARS-CoV-2 Nucleocapsid Protein is produced by E
<b>Conjugation</b>	Unconjugated
<b>Target</b>	SARS-CoV-2 Nucleocapsid
<b>Form/Appearance</b>	Lyophilized from a 0.2 µm filtered solution of PB, pH7.4
<b>Storage</b>	Ambient
<b>Tag</b>	C-His
<b>Note</b>	For research use only.
<b>Expression Region</b>	E.coli
<b>Uniprot ID</b>	<b>PODTC9</b>
<b>Expiration Date</b>	6 months from date of receipt.