

Product Datasheet

Normal Bovine Serum orb700576

Description	Normal Bovine Serum
Reactivity	Bovine
Hazard Information	Non-Toxic
Form/Appearance	Lyophilized
Storage	Store lyophilized material at 2-8°C. For long term storage after reconstitution, prepare small aliquots and store at -20°C . For storage at 2-8°C
Note	For research use only.
Purity	Lipid extracted and dialyzed against 10 mM Sodium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Expiration Date	12 months from date of receipt.