

## Product Datasheet

### Normal Bovine Serum orb700577

<b>Description</b>	Normal Bovine Serum
<b>Reactivity</b>	Bovine
<b>Hazard Information</b>	Non-Toxic
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	Store lyophilized material at 2-8°C. For long term storage after reconstitution, prepare small aliquots and store at -20°C . For storage at 2-8°C
<b>Note</b>	For research use only.
<b>Purity</b>	Lipid extracted and dialyzed against 10 mM Sodium Phosphate, 0.15 M Sodium Chloride, pH 7.2
<b>Expiration Date</b>	12 months from date of receipt.