

## Product Datasheet

### Collagenase/Elastase orb420143

|                           |  |
|---------------------------|--|
| <b>Description</b>        | Collagenase/Elastase   |
| <b>Hazard Information</b> | Non-Toxic  |
| <b>Target</b>             | Collagenase/Elastase   |
| <b>Form/Appearance</b>    | Lyophilized  |
| <b>Concentration</b>      | >22500 units collagenase/30 units elastase per vial  |
| <b>Storage</b>            | Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles. |
| <b>Note</b>               | For research use only.   |
| <b>Expiration Date</b>    | 12 months from date of receipt.  |