

Product Datasheet

Custom peptide synthesis (16 amino acids) orb98001

Description

Biorbyt offers a dedicated service for the production of high quality peptides. We work with an experienced team of organic chemists that work on **routine and complex synthesis**. We have **strict quality procedures** and provide our customers with the QC data relating to their products.

For simple synthetic requirements (6-30mer), place an order by adjusting the quantity (Qty default is 1) and email us your sequence with the order number. For more complex synthetic requirements, different purity levels and other related questions, please get in touch.

Price displayed is for each residue of a 6-30 mer peptide.

Hazard Information

Non-Toxic

Target

Custom peptide synthesis (16 amino acids)

Storage

Store at -20 °C. Avoid repeated freeze/thaw cycles.

Note

For research use only.

Purity

85% purified

Expiration Date

12 months from date of receipt.