

Product Datasheet

ApoA4 antibody orb704289

Description	Rabbit polyclonal antibody to ApoA4
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Custom conjugations	Unconjugated, Biotin, FITC, HRP, CF647, CF488A, CF594, CF405M
Tested Applications	ELISA, WB
Immunogen	KLH conjugated synthetic peptide derived from the region aa31-110 of human ApoA4
Target	ApoA4
Form/Appearance	Liquid: 0.01M TBS (pH7.4), 0.5% BSA, 0.015% ProClin™ 300 and 25% glycerol
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml - 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purified by Protein A
MW	41 kDa
Uniprot ID	P06727 (Human)
Entrez	337
Dilution Range	WB: 1:500-2000, ELISA: 1:500-1000
Expiration Date	12 months from date of receipt.