

Product Datasheet

Chlamydia trachomatis Virulence plasmid protein pGP3-D protein (orb704784)

Description	Recombinant Chlamydia trachomatis Virulence plasmid protein pGP3-D
Reactivity	Bacteria
Endotoxins	Not test.
Target	N/A
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 10xHis-tagged
Note	For research use only.
Protein Sequence	GNSGFYLYNTQNCVFADNIKVGQMTEPLKDQQIILGTTSTPVAAKMTASD GISLTVSNNPSTNASITIGLDAEKAYQLILEKLGDQILGGIADTIVDSTV QDILDKITTDPSLGLLKAFNNFPITNKIQCNGLFTPRNIETLLGGTEIGK FTVTPKSSGSMFLVSADIIASRMEGGVVLALVREGDSKPYAISYGYSSGV PNLCSLRTRIINTGLTPTTYSLRVGGLESGVWVWNALSNGNDILGITNTS NVSFLEVIPQTNA
Purity	Greater than 85% as determined by SDS-PAGE.
MW	30.3 kDa
Source	Yeast
Expression Region	2-264aa
Uniprot ID	POCE18
Expiration Date	6 months from date of receipt.