

Product Datasheet

Human ARF6 protein orb424150

Description	Recombinant of human ARF6 protein
Reactivity	Human
Conjugation	Unconjugated
Target	ARF6
Form/Appearance	Sterile Filtered clear solution. The ARF6 protein solution (1mg/1ml) is formulated in 20mM Tris-HCl buffer (pH 8.0), 2mM DTT, 0.2mM PMSF and 20% glycerol.
Storage	Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.
Note	For research use only.
Protein Sequence	MGSSHHHHHH SSGLVPRGSH MGKVLKIFG NKEMRILMLG LDAAGKTTIL YKLKLGQSVT TIPTVGFNVE TVTYKNVKFN VWDVGGQDKI RPLWRHYTG TQGLIFVVDC ADRDRIDEAR QELHRIINDR EMRDAILIF ANKQDLPDAM KPHEIQEKLK LTRIRDRNWY VQPSCATSGD GLYEGLTWLT SNYKS
Purity	>95% as determined by SDS-PAGE.
Source	Escherichia Coli.
Expiration Date	6 months from date of receipt.